PROVITEMP

TEMPORARY CEMENT

FR

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EN

Non-eugenol, resin based temporary cement containing Fluoride, Chlorhexidine and Potassium Nitrate.

INDICATIONS:

Temporary cementation of crowns, bridges, inlays, onlays and splints.

PROPERTIES:

ProviTemp has excellent flow to permit the restoration to be easily and completely seated. • It is strong enough to withstand the stresses of mastication, yet permits easy removal of the restoration. • ProviTemp does not contain essential oils so it will not harm temporary crowns or interfere with permanent cementation.

INSTRUCTIONS FOR USE:

- 1) Complete tooth preparation and fabrication of the temporary restoration before extruding the mixed cement.
- 2) Dry and isolate the prepared tooth. Be careful not to desiccate or over-dry the preparation.
- 3) Remove the cap from double barrel automix syringe and set aside. Extrude a very small portion of the unmixed material to insure unobstructed flow from both the base and catalyst components of the respective barrels.
- 4) Place the mixing tip onto syringe and lock in place with a quarter turn.
- 5) Extrude a thin layer of the automixed cement onto the cleaned inside surface of the temporary restoration. You have approximately 45 seconds of working time after extruding the mixed cement.
- 6) Firmly seat the temporary restoration in place. Initial gel-set occurs within 1.5 to 2 minutes after seating. Cleaning of excess cement should be initiated immediately after this gel-set phase.
- 7) Final cement set occurs between 3 and 4 minutes after seating of the temporary restoration.

TECHNICAL DATA:

Working time: 45 s*
Setting time: 3-4 min*
*These indications are obtained at the ideal mixing conditions of 24 °C (75 °F) and 40 % - 60 % relative humidity.

RECOMMENDATIONS/ WARNINGS:

Some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use and remove from patient. The patient should consult a doctor.

ESSENTIAL INGREDIENTS:

Fluoride • Chlorhexidine • Potassium Nitrate • Methacrylates • Urethan diméthacrylate • Polymerisation activator.

STORAGE:

Keep refrigerated. For dental use only.

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