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## ▶ High precision silicone-based addition-curing impression material

**Hydrospeed HD** is a complete range of high precision Polyvinylsiloxane addition-type silicone impression materials. Any of the products in the range may be combined with any other for use in one-step technique or in two-step technique (putty-wash technique).

## ▶ Properties

The **Hydrospeed HD** range is 100% biocompatible and nontoxic; it offers a variety of hardnesses and viscosities, and so covers all impression-taking techniques. Its hydro-compatibility ensures perfect reproduction of details. Being thixotropic, the material gives optimal flow in the sub-gingival zones. The impressions made do not go out of shape; they are highly accurate due to the excellent mechanical properties (elasticity, dimensional stability, &c.) of the silicone elastomers used. The materials of the **Hydrospeed HD** range also have good resistance to light and chemicals.

## ▶ Indications

The **Hydrospeed HD** range may be used for all type of impression where accuracy is needed: crowns • bridges • inlays • onlays • veneers • splints • orthodontic impressions • antagonists.

## ▶ The Hydrospeed HD range

| Name of Product                       | Classification under ISO 4823 | Setting time | Colour    | Packaging                | Accessories                  |
|---------------------------------------|-------------------------------|--------------|-----------|--------------------------|------------------------------|
| Hydrospeed HD Putty Hard              | Type 0 - Putty                | Normal       | Dark pink | 2 x 300 ml jars          | 2 spoons                     |
| Hydrospeed HD Putty Hard Quick        | Type 0 - Putty                | Quick        | Green     | 2 x 300 ml jars          | 2 spoons                     |
| Hydrospeed HD Magnum Putty Hard       | Type 0 - Putty                | Normal       | Blue      | 1 x 380 ml cartridge 5:1 | --                           |
| Hydrospeed HD Magnum Putty Hard Quick | Type 0 - Putty                | Quick        | Green     | 1 x 380 ml cartridge 5:1 | --                           |
| Hydrospeed HD Putty Soft              | Type 0 - Putty                | Normal       | Yellow    | 2 x 300 ml jars          | 2 spoons                     |
| Hydrospeed HD Putty Soft Quick        | Type 0 - Putty                | Quick        | Sky blue  | 2 x 300 ml jars          | 2 spoons                     |
| Hydrospeed HD Magnum Body Heavy       | Type 1 - Heavy                | Normal       | Yellow    | 1 x 380 ml cartridge 5:1 | --                           |
| Hydrospeed HD Magnum Body Heavy Quick | Type 1 - Heavy                | Quick        | Sky blue  | 1 x 380 ml cartridge 5:1 | --                           |
| Hydrospeed HD Regular Body            | Type 2 - Regular              | Normal       | Pink      | 2 x 50 ml cartridges 1:1 | 12 yellow EcoMix mixing tips |
| Hydrospeed HD Regular Body Quick      | Type 2 - Regular              | Quick        | Grey      | 2 x 50 ml cartridges 1:1 | 12 yellow EcoMix mixing tips |
| Hydrospeed HD Light Body              | Type 3 - Light                | Normal       | Green     | 2 x 50 ml cartridges 1:1 | 12 yellow EcoMix mixing tips |
| Hydrospeed HD Light Body Quick        | Type 3 - Light                | Quick        | Orange    | 2 x 50 ml cartridges 1:1 | 12 yellow EcoMix mixing tips |

| Name of Product                       | Packaging                | Setting time | Colour    | Hardness (Shore A) | Working time <sup>1</sup> | Setting time in the mouth <sup>2</sup> | Total time <sup>3</sup> |
|---------------------------------------|--------------------------|--------------|-----------|--------------------|---------------------------|--|-------------------------|
| Hydrospeed HD Putty Hard              | 2 x 300 ml jars          | Normal       | Dark pink | 70 A               | 2 min 20                  | 3 min 30                               | 5 min 50                |
| Hydrospeed HD Putty Hard Quick        | 2 x 300 ml jars          | Quick        | Green     | 70 A               | 1 min 30                  | 2 min 35                               | 4 min 05                |
| Hydrospeed HD Magnum Putty Hard       | 1 x 380 ml cartridge 5:1 | Normal       | Blue      | 68 A               | 2 min 20                  | 3 min 30                               | 5 min 30                |
| Hydrospeed HD Magnum Putty Hard Quick | 1 x 380 ml cartridge 5:1 | Quick        | Green     | 68 A               | 1 min 30                  | 2 min 35                               | 4 min 05                |
| Hydrospeed HD Putty Soft              | 2 x 300 ml jars          | Normal       | Yellow    | 65 A               | 2 min 20                  | 3 min 30                               | 5 min 50                |
| Hydrospeed HD Putty Soft Quick        | 2 x 300 ml jars          | Quick        | Sky blue  | 65 A               | 1 min 30                  | 2 min                                  | 3 min 30                |

### ▶ PUTTY

1 - Working time: typical duration of mixing and handling before placing in the mouth.  
 2 - Setting time in the mouth: typical time for optimum material curing at 32°C. When the temperature is above 23°C the setting time will be shorter, and conversely.  
 3 - Total time: typical time from the beginning of mixing until complete material cure.

### Dosage and mixing

Take equal quantities of the two component materials, using the right spoon for each. To reduce any risk of polymerising the materials prematurely, the colour of each spoon matches that of the appropriate material. Mix with the fingers for approximately 30 secs until the colour is uniform throughout. Close the jars immediately after use. Be sure to correctly replace the lids and the spoons used to measure doses of base material and catalyst.


## ▶ BODY HEAVY TRAY AND WASH MATERIAL

| Name of Product                       | Packaging                | Setting time | Colour   | Hardness (Shore A) | Viscosity | Working time <sup>1</sup> | Setting time in the mouth <sup>2</sup> | Total time <sup>3</sup> |
|---------------------------------------|--------------------------|--------------|----------|--------------------|-----------|---------------------------|--|-------------------------|
| Hydrospeed HD Magnum Body Heavy       | 1 x 380 ml cartridge 5:1 | Normal       | Yellow   | 62 A               | Heavy     | 2 min 20                  | 3 min 30                               | 5 min 50                |
| Hydrospeed HD Magnum Body Heavy Quick | 1 x 380 ml cartridge 5:1 | Rapid        | Sky blue | 62 A               | Heavy     | 1 min 30                  | 2 min 35                               | 4 min 05                |
| Hydrospeed HD Regular Body            | 2 x 50 ml cartridges 1:1 | Normal       | Pink     | 45 A               | Regular   | 2 min 30                  | 3 min 30                               | 6 min                   |
| Hydrospeed HD Regular Body Quick      | 2 x 50 ml cartridges 1:1 | Rapid        | Grey     | 45 A               | Regular   | 1 min 30                  | 2 min 35                               | 4 min 05                |
| Hydrospeed HD Light Body              | 2 x 50 ml cartridges 1:1 | Normal       | Green    | 53 A               | Light     | 2 min 30                  | 3 min 30                               | 6 min                   |
| Hydrospeed HD Light Body Quick        | 2 x 50 ml cartridges 1:1 | Rapid        | Orange   | 53 A               | Light     | 1 min 30                  | 2 min 35                               | 4 min 05                |

1 - Working time: typical duration of mixing and handling before placing in the mouth.  
 2 - Setting time in the mouth: typical time for optimum material curing at 32°C. When the temperature is above 23°C the setting time will be shorter, and conversely.  
 3 - Total time: typical time from the beginning of mixing until complete material cure.

## ▶ Dispensing gun

**Hydrospeed HD** cartridges are compatible with all mixing dispensers on the market, including the ITENA 1:1 mixing dispenser. Standard market mixing tips (4.2mm diameter) are compatible with the **Hydrospeed HD** cartridges.

1 Insert the cartridge in the blue guide ring.

2 Lock in place with an anticlockwise twist.

3 Remove the sealing cap and extrude a little product until the two openings are producing material at the same rate.

4 Fit an intra-oral tip (if necessary) onto the yellow EcoMix tip and fix the assembly to the cartridge with a 90-degree clockwise twist.


5 As the trigger is squeezed the material is extruded and automatically mixed in the right proportion; it may be applied directly in the mouth to the contact point of the preparation(s). Maintain steady, uninterrupted pressure. Flow ceases as soon as the pressure is relaxed.

6 Remove the intra-oral tip. Disinfect the dispensing gun, cartridge and tip before putting them away.

Do not remove the yellow EcoMix tip from the cartridge after use; leave it in place until next time, when it will be replaced. When changing the yellow EcoMix tip check that the cartridge openings are not blocked so that the material will be able to flow freely.

## ▶ Cleaning and disassembly

The mixing dispenser and **Hydrospeed HD** cartridges can be cleaned using antiseptic solution or wipes.

1 Turn the blue guide ring all the way.

2 Keep the dispenser clamp uppermost and remove the piston.

3 The dispenser is disassembled.

## ▶ Sterilisation

The Dispenser handle and piston can be autoclaved up to a temperature of 135°C/275°F.

## The Magnum CARTRIDGE

### ▶ Description of the product

The cartridge delivers base material and catalyst at a ratio of 5:1. **Hydrospeed HD Magnum** cartridges are wholly compatible with mixing machines meeting the ISO 4823 standards such as: MixStar-eMotion; DMG • PowerMix; Kaniedenta • Sympress; Kettenbach • Plug & Press; Kettenbach • Duomix II; Dentsply • Modulmix; Zhermack • Dynamix speed; Heraeus Kulzer • Sympress; Renfert • Pentamix 2; 3M • Pentamix 3; 3M \*

\*Hydrospeed HD Magnum cartridges only work at low speed in the 3M Pentamix 3 mixing machine. This in no way affects the quality of the mixture.

**Hydrospeed HD Magnum** cartridges are compatible with dynamic mixing tips and their associated fixation rings (including ITENA tips), except for 3M tips specifically designed for 3M Pentamix cartridges.

### ▶ Loading the Magnum cartridge

See the manufacturer's instructions for your mixing machine.

1. Hold the **Hydrospeed HD Magnum** cartridge upright, extrusion end uppermost.
2. Take hold of the tab of the protective plastic end cap and pull gently upwards.
3. Remove the cap, holding the **Hydrospeed HD Magnum** cartridge firmly in the hand.
4. Insert the cartridge in the mixing machine, following the manufacturer's instructions.
5. Extrude a little material to start. Repeat this operation at the start of every use. Wipe the end of the cartridge to prevent any contamination between base material and catalyst. On new **Hydrospeed HD Magnum** cartridges the first centimetres of material extruded must be discarded.
6. Fit a suitable mixing tip to the cartridge.
7. Make sure that the hexagonal sleeve is correctly engaged at the start.
8. Put a suitable fixation ring in position for the tip. Give a quarter turn to lock.
9. Extrude the silicone into a suitable impression tray.
10. After use disinfect following the manufacturer's recommendations; keep the used mixing tip in place on the cartridge until the next use.

### ▶ Cleaning the Hydrospeed HD Magnum cartridge

**Hydrospeed HD Magnum** cartridges can be cleaned using solution or disinfectant wipes.

## ▶ IMPRESSION TECHNIQUES

### ➔ One-step technique

1. Put the Fluid impression material around the preparation directly in the mouth.
  - 1.1 An intra-oral syringe may also be used, the material being extruded directly into the syringe.
2. Fill the impression tray with the Putty impression material.
3. Put the impression tray into the mouth; leave for the chosen material's setting time before removing.
4. When the material is polymerised open the patient's mouth and gently remove it. Remove any excess material from undercut areas.

### ➔ Two-steps technique (Wash technique)

1. Fill the impression tray with the Putty impression material.
2. Put the impression tray into the mouth and remove after the paste has set.
3. Clean and prepare the primary impression obtained.
4. Use a water jet to clean the area to be registered; dry carefully.
5. Apply the Fluid impression material to the area to be registered directly in the mouth.
  - 5.1 An intra-oral syringe may also be used, the material being extruded directly into the syringe.
6. Squirt Fluid impression material onto the occlusal surfaces of the adjacent teeth as well, and/or in the impression tray containing the primary impression.
7. Re-insert the impression tray in the mouth, and keep it in position for the chosen material's setting time.
8. When the material is polymerised open the patient's mouth and gently remove it. Remove any excess material from undercut areas.

### ▶ Impression tray

All impression tray used for precision impressions are suitable for **Hydrospeed HD** products.

### ▶ Disinfecting the impression

The impression may be cleaned using impression-disinfecting solution or wipes. Follow the manufacturer's instructions, particularly as regards time soaking in the solution so as not to risk degrading the impression. Then rinse the impression under running water to remove any trace of disinfecting solution.

### ▶ Model preparation

All dental plaster on the market can be used for pouring the impression. Do not pour less than 30 minutes after removal from the mouth. Long-term stability: the impression is guaranteed to last for one month or more under proper storage conditions. Dimensional changes are no greater than 0.2% (ISO 4823).

### ▶ Galvanisation

Impressions made using the **Hydrospeed HD** can be silver or copper-plated.

### ▶ Precautions for use

Never mix Polyvinylsiloxane addition-type materials with silicones by condensation or with polyether impression materials. Any such combination will delay or prevent setting. Use vinyl or nitrile gloves in preference to powdered latex, which hinders the setting of silicone-based impression materials. Do not allow impression materials to come into contact with solvents.

### ▶ Contraindications

Use of this product can cause an allergic reaction in rare cases. Discontinue the procedure immediately and consult a doctor.

### ▶ Storage and stability

**Hydrospeed HD** impression materials must not be used beyond their stated expiration date. They must be kept in a dry place and store at -5°C - 30°C.

**For dental use only. Keep out of the reach of children**

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